

Hingmed 24-Hour Ambulatory Blood Pressure Monitor

Component option:

Item	Cuff Option
1	Small adult cuff (18-26cm)
2	Middle adult cuff (22-32cm)
3	Large adult cuff (22-36cm)
4	Extra large adult cuff (33-43cm)



WBP-02A includes:

Parts	Number
Hingmed ABPM software U-disk	1
Hingmed ABPM USB data cable	1
Quick battery charger	1
Large adult cuff (22-36 cm)	1
Silicone cover	1
Disposable sweat absorbable inner cuff	1
Hingmed ABPM bag	1

WBP-02A Product Parameters

Name	24-Hour ABPM
Model	WBP-02A
Technology	Oscillometric
Systolic	40-260 mmHg
Diastolic	20-210 mmHg
Heart Rate	40-200 BPM
Accuracy	± 3 mmHg
Pressure Transducer	America Freescale
Battery	1000mah Li battery,support 200 times measurement
Memory	Saving 300 readings Maximum
Calibration	At least once per two years
Safety	Max inflation Pressure290mmHg; Valve open for deflation automatically upon power shutting off; Max measurement circle 120s.
Intercals Option	More individual programming intervals options (5,10,15,20,30,45,60,90、120min)
Size	119 x 52 x 21 mm
Weight	105g, including battery
Working Condition	Tem: 5°C ~ 40°C; RH: 10%~95%; Atmospheric pressure: 70Kpa~106Kpa;
Storage Condition	Tem:-20°C~+55°C; RH less than 95%; Atmospheric pressure: 70Kpa~106Kpa;
Data Connection	USB data cable

HINGMED

WBP-02A

24-Hour Ambulatory Blood Pressure Monitor



Shenzhen Hingmed Medical Instrument Co.,Ltd

Address: 4F,Zhonghangfeixiang Building,No.371,Guangshen Road,Bao'an District,Shenzhen,China.  
Mobile/WhatsApp/Wechat: +86 -18681556084  
Website:http://www.hingmedical.com  
http://en.hingmed.com  
Email:sales@hingmed.com



Hingmed Medical focuses on blood pressure and respiratory monitoring technology and products. The intelligent pulse wave inflation measurement technology provides a very comfortable BP measurement experience. WBP-02A is a wearable 24-Hour ABPM without air tube, provides very comfortable 24-Hour BP monitoring experience. WBP-02A is ESH IP 2010 validated, recommended by ESH and BIHS, listed in Dabl Educational Trust and Stride BP.

# HINGMED

- Small, Light, Wearable.
- Clinical level, accurate reading.
- Motion tolerance technology, good measurement consistency.
- Body position record, good help for doctor's analysis.



- Double pressure transducer, double protection.
- BP limits alarm.
- Low power consumption, one fully charging support 200 time measurement.



- Microsoft Windows 2003、Windows XP or higher version application environment.
- Support maximum 72 hours data tendency chart.
- Analysis of morning surge, smoothness index, and AASI.
- White coat analysis, remarkable color for white coat period.
- Various BP reading analysis, statistics.
- BP comparison of same patient's different measurement.
- Variation correlation analysis of DIA vs SYS., HR vs SYS.
- PDF format report.
- User defined report content, printing review.

Hingmed WBP-02A ABPM System

