

## COTTON HYDROPHILIC GAUZE - STERILE

### STERILE COTTON HYDROPHILIC GAUZE

Folded hydrophilic gauze 100% pure cotton, suitable for sterile medical dressings, to buffer and protect wounds. It allows skin perspiration. Individually packed and latex free.

Instructions and box in GB, FR, IT, ES, DE, GR

MADE IN ITALY



| GIMA code | Threads | COTTON HYDROPHILIC GAUZE - STERILE | Minimum order              |
|-----------|---------|------------------------------------|----------------------------|
| 34768     | 12/8    | Size 18x40 cm                      | carton of 48 box of 12 pcs |
| 34769     | 12/8    | Size 36x40 cm                      | carton of 48 box of 12 pcs |

MADE IN ITALY



Instructions and box in GB, IT

| GIMA code | Threads | COTTON HYDROPHILIC GAUZE - STERILE | Minimum order                |
|-----------|---------|------------------------------------|------------------------------|
| 34765     | 12/8    | Size 10x10 cm                      | carton of 100 bags of 25 pcs |
| 34766     | 12/8    | Size 15x15 cm                      | carton of 100 bags of 25 pcs |
| 34767     | 12/8    | Size 20x20 cm                      | carton of 100 bags of 25 pcs |

## LAPAROTOMIC SPONGE GAUZE



## NON WOVEN FABRIC GAUZE - STERILE

### NON WOVEN FABRIC GAUZE WITHOUT X-RAY

Non invasive device to be used as a mechanical barrier, for compression or for absorption of exudates. Intended to manage the micro environment of a wound. Maximum absorbance, it does not adhere to the wound. Made of viscose and polyester. Medical device class risk I

MADE IN EUROPE

Individually packed.



| GIMA code | NON WOVEN FABRIC GAUZE - STERILE | Minimum order              |
|-----------|----------------------------------|----------------------------|
| 35026     | Size 18x40 cm                    | carton of 50 box of 12 pcs |
| 35027     | Size 36x40 cm                    | carton of 40 box of 12 pcs |

MADE IN EUROPE



| GIMA code | NON WOVEN FABRIC GAUZE - STERILE | Minimum order            |
|-----------|----------------------------------|--------------------------|
| 35024     | Size 10x10 cm                    | carton of 60 box of 100* |

\*4 bags of 25 pcs

### LAPAROTOMIC SPONGE GAUZE - not sterile

100% cotton gauze made of 4 layers (8 layers for 59394) of 20 threads/cm², sewn on sides and highly absorbent. Barium X-ray thread woven into the gauze, made of biocompatible polypropylene, with 60% barium sulphate. Class IIa Medical Device according to 93/42/CEE, 2007/47/CEE Latex free.

| GIMA code | Plies | Threads | LAPAROTOMIC SPONGE GAUZE | Minimum order |
|-----------|-------|---------|--------------------------|---------------|
| 59390     | 4     | 20      | Size 10x60 cm            | box of 600    |
| 59391     | 4     | 20      | Size 20x20 cm            | box of 600    |
| 59392     | 4     | 20      | Size 30x30 cm            | box of 600    |
| 59393     | 4     | 20      | Size 40x40 cm            | box of 250    |
| 59394     | 8     | 20      | Size 40x40 cm            | box of 250    |
| 59395     | 4     | 20      | Size 50x50 cm            | box of 250    |

## 100% COTTON GAUZE ROLLS




### GAUZE ROLLS

100% cotton gauze rolls with 20 threads/sq cm according to European Pharmacopoeia and EN 14079. X-Ray detectable thread (approx. 60% barium sulphate) 80 cm width, folded 8 to 10 cm ply. Multilingual box: GB, FR, IT, ES, PT, DE, NL, SE.

| GIMA code | Extended width | Folded width | Plies | Threads/sq cm |
|-----------|----------------|--------------|-------|---------------|
| 59410     | 80 cm x 1 m    | 10 cm x 1 m  | 8     | 20            |
| 59411     | 80 cm x 2 m    | 10 cm x 2 m  | 8     | 20            |
| 59412     | 80 cm x 10 m   | 10 cm x 10 m | 8     | 20            |

## 100% COTTON RIBBON GAUZE

- Latex and PVC free
- Absorption of fluids inside cavities
- Wound care during surgery



Multilingual box: GB, FR, IT, ES, PT, DE, NL, SE

| GIMA code | GIMA code | 100% COTTON RIBBON GAUZE weight 33 g/sq. m 22 threads/sq. m | Minimum order |
|-----------|-----------|-------------------------------------------------------------|---------------|
| -         | X-Ray*    |                                                             |               |
| 59420     | 59430     | 5 m x 1.25 cm                                               | box of 24     |
| 59421     | 59431     | 5 m x 2.5 cm                                                | box of 12     |
| 59422     | 59432     | 5 m x 5 cm                                                  | box of 12     |
| 59423     | 59433     | 5 m x 7.5 cm                                                | box of 8      |
| 59424     | 59434     | 5 m x 10 cm                                                 | box of 6      |
| 59440     | -         | 5 m x 1.25 cm - sterile                                     | box of 60     |
| 59441     | -         | 5 m x 2.5 cm - sterile                                      | box of 30     |
| 59442     | -         | 5 m x 5 cm - sterile                                        | box of 26     |

\*X-Ray: approx. 60% barium sulphate contents