


Certificate

Polymer Technology Systems Inc. (PTS Diagnostics)

using

A1CNow+, A1CNow SELF CHECK

participated in the IFCC HbA1c Certification Programme to demonstrate traceability to the IFCC Reference Measurement Procedure and performed as shown below.



Total Error =	3.3	mmol/mol
Bias =	-0.8	mmol/mol
Imprecision =	2.5	%
Grade =	Bronze	

Criteria derived from the IFCC model for Quality Targets HbA1c (Clin Chem 2015;61 : 752-59)

Date of Certification : 01 January 2025

Date of Expiry : 01 January 2026

G. Siebelder
IFCC Network Coordinator
G. Siebelder