



Otoscope EOC 700

— Full Color Display

Full HD Image



Ð

Tympanic Membrane & Throat Inspection

Dual inspection modes of tympanic membrane and throat provide the best diagnostic experience for the doctors.



















Emergency



Ward round

Telemedicine

Education

Veterinary

Model	horus SCOPE _® – Otoscope EOC 700
Focus Type	Manual focus
Image Resolution	1920 x 1080 pixels
LCD Monitor	3.5" TFT-LCD
Image Format	JPEG (photograph) and H.264 (video)
Interface	Mini USB
File Transfer	Mini USB port to PC
File Storage	Memory card, supports 8GB to 32GB by FAT32 format
Power Source	Rechargeable lithium battery 3.6V/3350mAh The battery is in compliance with IEC 62133 2:2017
External Power	Source: 100~240 VAC, 50/60 Hz
Power Adapter Spec	Input spec.: 100~240 VAC, 50/60Hz, 0.2A Output spec.: 5V DC, 1.2A
Operating Time	3 hours
Charging Time	5 hours



Medimaging Integrated Solution Inc.

+886 3 579-8860 +886 3 579-8011 6F, No. 11, Park Avenue II, Hsinchu Science Park, Hsinchu 300091, Taiwan, R.O.C. www.miis.com.tw sales@miis.com.tw

