

Universal Wipes

Efficacy Data



Bacteria example	Test	Condition	Contact time
<i>Enterococcus hirae</i>	EN 16615	Dirty	30 seconds
	EN 13727	Dirty	10 seconds
<i>Pseudomonas aeruginosa</i>	EN 16615	Dirty	30 seconds
	EN 13727	Dirty	10 seconds
<i>Staphylococcus aureus</i>	EN 16615	Dirty	30 seconds
	EN 13727	Dirty	10 seconds
<i>Acinetobacter baumannii</i> (Carbapenem resistant)	EN 13727	Dirty	10 seconds
<i>Enterococcus faecium</i> (Vancomycin resistant - VRE)	EN 13727	Dirty	10 seconds
<i>Escherichia coli</i> k12	EN 13727	Dirty	10 seconds
<i>Klebsiella pneumoniae</i> (Carbapenem resistant)	EN 13727	Dirty	10 seconds
<i>Leptospira interrogans</i>	EN 13727	Dirty	30 seconds
Methicillin-resistant <i>Staphylococcus aureus</i> (MRSA)	EN 13727	Dirty	10 seconds
<i>Stenotrophomonas maltophilia</i>	EN 13727	Dirty	30 seconds

Yeast example	Test	Condition	Contact time
<i>Candida albicans</i>	EN 16615	Dirty	30 seconds
	EN 13624	Dirty	10 seconds
<i>Candidozyma auris</i>	EN 16615	Dirty	30 seconds

Virus example	Test	Condition	Contact time
Adenovirus	EN 16777	Dirty	30 seconds
	EN 14476	Dirty	30 seconds
Norovirus	EN 16777	Dirty	30 seconds
	EN 14476	Dirty	30 seconds
Coronavirus	EN14476	Dirty	30 seconds
Herpes simplex virus	EN 14476	Dirty	30 seconds
Hepatitis B	EN 14476	Dirty	30 seconds
Hepatitis C	EN 14476	Dirty	30 seconds
HIV	EN 14476	Dirty	30 seconds
Influenza	EN 14476	Dirty	30 seconds
Measles virus	EN 14476	Dirty	30 seconds
Mumps virus	EN 14476	Dirty	30 seconds
Papillomavirus (HPV)	EN 14476	Dirty	30 seconds
Papovavirus (SV40)	EN 14476	Dirty	30 seconds
Rotavirus	EN 14476	Dirty	30 seconds
Vaccinia virus	EN 16777	Dirty	15 seconds
	EN 14476	Dirty	15 seconds

All efficacy data is validated in independent UKAS-accredited, ISO-certified laboratories.